
PCN Number: 08/2024

Title: **Changed cable for PowerPack for Microscopy with pulse camera**

Product Family: Powerpack for Microscopy with pulse

Release date: October 10, 2024

Contact: Basler Support

Type of Change:

Miscellaneous Firmware PC Software

Product design Process

Reason for Change

Due to the availability of the USB 3.0 cable which is currently included in the microscopy PowerPacks with pulse cameras, we are forced to change the cable to another one.

Change Description

The previous cable was **angled backwards** on the camera side to have the cable directly at the backside of the microscope. (see picture)



Figure 1: previous USB Cable

To keep the deliverability on these products we decided to replace the previous cable with a **straight** USB 3.0 cable from our existing portfolio.

Basler AG
Europe, Middle East, Africa
Tel. +49 4102 463 515
Fax +49 4102 463 599
support.europe@baslerweb.com
www.baslerweb.com

The Americas
Tel. +1 610 280 0171
Fax +1 610 280 7608
support.usa@baslerweb.com

Asia-Pacific Region
Tel. +65 6367 1355
Fax +65 6367 1255
support.asia@baslerweb.com

This may cause a slightly different handling of the cable. In all other concerns the functionality is fully given.

The new cable has screw locks that the previous cable didn't have. However, the pulse camera doesn't support screw locks, so they are not usable and can be ignored.

Affected models and expected first production / shipping date

Model	Revision	Production /shipping Date
<i>PowerPack for Microscopy with Microscopy pulse 1.2 MP</i> PowerPack puA1280-54ucMIC - 107226	09	CW42 / 2024
<i>PowerPack for Microscopy with Microscopy pulse 2.0 MP</i> PowerPack puA1920-30ucMIC - 107227	09	CW42 / 2024
<i>PowerPack for Microscopy with Microscopy pulse 5.0 MP</i> PowerPack puA2500-14ucMIC - 107228	09	CW42 / 2024
<i>PowerPack for Microscopy with Microscopy pulse 3.3 MP</i> PowerPack puA2100-20ucMIC - 107245	09	CW42 / 2024

Additional Information

Basler may make changes to specifications and product descriptions at any time, without notice. The information in this document is subject to change without notice.

There is no claim to exchange or update the stock or the delivered goods due to the PCN. This applies to customer inventories as well as Basler's own inventory. Also there is no claim to receive the old revision after change.

Basler AG
Europe, Middle East, Africa
Tel. +49 4102 463 515
Fax +49 4102 463 599
support.europe@baslerweb.com
www.baslerweb.com

The Americas
Tel. +1 610 280 0171
Fax +1 610 280 7608
support.usa@baslerweb.com

Asia-Pacific Region
Tel. +65 6367 1355
Fax +65 6367 1255
support.asia@baslerweb.com